

• If a user is adding a new location to the database then they shall click on the Location button on the top menu bar.

Home Contents Lookup Add an Aliquot CoA Test Page Shipping Returns Aliquot Search Form Add a Box Upload CSV? Help!? Location

 The submenu for Locations will expand and the user will click on Add a New Location.




- Once the user clicks on the Add a New Location link they will be redirected to the Add a Location page.
- All fields are required when adding a new location.




• Once the user has filled out the four fields they can press the submit button and the new location will be added.

• If a user is wanting to check which locations are already in the database, which is a good idea before adding a new location. The user will click on **Locations List** in the submenu of Location.



 To search the page the user will be able to filter the page to narrow the displayed sites.



• If the user wants to download the site information to a csv document then they shall click on the **Export Table to CSV** button in the top left of the table.

