

- Contents Look Up page
- Use search if you know what contents you are looking for
- Use View All to scroll all contents.

Contents Lookup						
<div> Contents Lookup Freezer Inventory Add an Aliquot Aliquot Table Add a Box Upload CSV? </div>						
<div> <div>Search...</div> <div> Search Export Table to CSV View All Contents </div> </div>						
Box Type	Box ID	Contents	Insert Date	Inventory Status	Freezer	
10x10	001A	ANA Positives	April 23, 2015, 5 p.m.	1	240B	
9x9	001B	ANA Positives Negatives	April 23, 2015, 5 p.m.	1	240B	
10x10	002A	SLE Positives - QCP inc Positives	April 23, 2015, 5 p.m.	1	240B	
9x9	002B	SLE Positives	April 23, 2015, 5 p.m.	1	240B	
10x10	003A	Anti-CMV Pos+ IgG	April 23, 2015, 5 p.m.	1	240B	
9x9	004A	Anti-CMV Positive IgM Positives	April 23, 2015, 5 p.m.	1	240B	
10x10	005A	Anti-EBV Positive	April 23, 2015, 5 p.m.	1	240B	

- To search for specific contents either use the search bar or the dropdowns to filter your view.

Contents Lookup

001A

Search

Export Table to CSV

View All Contents

Box Type

Box ID

Contents

Insert Date

Inventory Status

Freezer

Box Type	Box ID	Contents	Insert Date	Inventory Status	Freezer
10x10	001A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
10x10	053A	HIV-1 Seroconversion Panels, Other, HIVSCP-001a,SV	April 23, 2015, 5 p.m.	0	Unknown

Contents Lookup

Search....

Search

Export Table to CSV

View All Contents

Box Type

Box ID

Contents

Insert Date

Inventory Status

Freezer

Box Type	Box ID	Contents	Insert Date	Inventory Status	Freezer
10x10	001A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
9x9	001B	ANA Positives	April 23, 2015, 5 p.m.	1	240B
10x10	002A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
9x9	002B	ANA Positives	April 23, 2015, 5 p.m.	1	240B
10x10	003A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
9x9	004A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
10x10	005A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
9x9	005B	ANA Positives	April 23, 2015, 5 p.m.	1	240B

- To export an excel file with the view the user has filtered, click on the 'Export Table to CSV'.

Ana positives

Search **Export Table to CSV** View All Contents

Box Type	Box ID	Contents	Insert Date	Inventory Status	Freezer
10x10	001A	ANA Positives	April 23, 2015, 5 p.m.	1	240B
9x9	001B	ANA Positives Negatives	April 23, 2015, 5 p.m.	1	240B

- Example of the exported content.

A	B	C	D	E	F
Box Type	Box ID	Contents	Insert Date	Inventory Status	Freezer
10x10	001A	ANA Positives	April 23: 2015: 5 p.m.	1	240B
9x9	001B	ANA Positives Negatives	April 23: 2015: 5 p.m.	1	240B

- Once a user has found the contents they are looking for they can click on any of the cells in the view and it will redirect them to the Aliquot and Box detail view.
- This is the main transaction page for the repository application.
- You can update the Box information, Export the Contents of the box to excel, checkout aliquots, and update aliquot information.

Aliquot and Box Detail											
<div>Export Table to CSV</div> <div>You are in box: 001A. If you would like to update the box 001A click--> Update</div>											
Select Aliquot	AQID	Cage	Cane	Stack	Freezer	Original ID	Grid	Label Information	Volume (mL)	Insert Date	Checked Out?
<input type="checkbox"/>	001A-A01- Update	1	01	A01	240B	00047303	A01	00047303	1.000	April 23, 2015, 5 p.m.	No
<input type="checkbox"/>	001A-A02- Update	1	01	A01	240B	00047303	A02	00047303	1.000	April 23, 2015, 5 p.m.	No
<input type="checkbox"/>	001A-A03- Update	1	01	A01	240B	00047303	A03	00047303	1.000	April 23, 2015, 5 p.m.	No
<input type="checkbox"/>	001A-A05- Update	1	01	A01	240B	00047281	A05	00047281	1.000	April 23, 2015, 5 p.m.	No

- To checkout an aliquot the user will click the check boxes of all the aliquots they require and click the 'check out' button

USER: datamgt
[Logout](#)

Aliquot and Box Detail

[Export Table to CSV](#)

You are in box: **001A**. If you would like to update the box **001A** click--> [Update](#)

Select Aliquot	AQID	Cage	Cane	Stack	Freezer	Original ID	Grid	Label Information	Volume (mL)	Insert Date	Checked Out?
<input checked="" type="checkbox"/>	001A-A01-Update	1	01	A01	240B	00047303	A01	00047303	1.000	April 23, 2015, 5 p.m.	No

- Once the button is pressed the page will update with a success message in green and the column 'Select Aliquot' will no longer show check boxes where the aliquot is checked out. Also the 'Checked Out?' column will change from No to Yes.

USER: datamgt
[Logout](#)

Aliquots [001A-A01] Checked Out Successfully!

Aliquot and Box Detail

[Export Table to CSV](#)

You are in box: **001A**. If you would like to update the box **001A** click--> [Update](#)

Select Aliquot	AQID	Cage	Cane	Stack	Freezer	Original ID	Grid	Label Information	Volume (mL)	Insert Date	Checked Out?
	001A-A01-Update	1	01	A01	240B	00047303	A01	00047303	1.000	April 23, 2015, 5 p.m.	Yes