

- If a user needs to dispose of a sample then while in the Aliquot and Box Detail view they should press the dispose button. This will redirect the user to the Aliquot Dispose view.

Aliquot and Box Detail

[Export Table to CSV](#)

You are in box: **001A**. If you would like to update the box **001A** click--> [Update](#)

[Click here to Dispose of Aliquots - \[Dispose Aliquots\]\(#\)](#)

- Once the user is in the Aliquot Dispose view then they can select the aliquots they want to dispose of and press the disposal button on the left hand side of the table

Dispose	<input checked="" type="checkbox"/>	001A-testingagain123- Update	1.500	No
	<input checked="" type="checkbox"/>	001A-testingMulti-0- Update	None	No
	<input checked="" type="checkbox"/>	001A-testingmulti-11- Update	1.000	No
	<input type="checkbox"/>	001A-testingmulti-12- Update	1.000	No
	<input type="checkbox"/>	001A-testingmulti-13- Update	1.000	No
	<input type="checkbox"/>	001A-testingmulti-14- Update	1.000	No
	<input type="checkbox"/>	001A-testingmulti-15- Update	1.000	No

- After the dispose button is pressed a success message at the top will display all of the aliquots that were disposed.

Aliquots ['001A-testingagain123', '001A-testingMulti-0', '001A-testingmulti-11'] Disposed Successfully!

Aliquot Disposal

[Export Table to CSV](#)

You are in box: **001A**. If you would like to update the box **001A** click--> [Update](#)

Dispose Aliquot	AQID	Volume (mL)	Already Disposed?
Dispose Aliquot	AQID	Volume (mL)	Already Disposed?
<input type="checkbox"/>	001A-A01- Update	1.000	No

- **Note:** In the dispose view the right hand most column will read 'yes' if the aliquot has been disposed of.
- **Note:** This view is also a running list of all of the samples that have been assigned to the box. It will note which ones were disposed or which ones are still active in the box. Aliquots that are still active the right hand most column will read 'no'.
- **Note:** Once the aliquots have been disposed of they will not be visible in any view other than the Aliquot Dispose view.

<input type="checkbox"/>	001A-testingagain123- Update	1.500	Yes
<input type="checkbox"/>	001A-testingMulti-0- Update	None	Yes
<input type="checkbox"/>	001A-testingmulti-11- Update	1.000	Yes
<input type="checkbox"/>	001A-testingmulti-12- Update	1.000	No