

- To update the aliquot the user will click on the 'Update' link next to it's corresponding Aliquot ID in the AQID column.
- This procedure is also how a user will check an aliquot back into the box.
- The process is sped up by knowing the Aliquot ID

Aliquot and Box Detail											
<a href="#">Export Table to CSV</a>											
You are in box: <b>001A</b> . If you would like to update the box <b>001A</b> click--> <a href="#">Update</a>											
Select Aliquot ▼	AQID ▼	Cage ▼	Cane ▼	Stack ▼	Freezer ▼	Original ID ▼	Grid ▼	Label Information ▼	Volume (mL) ▼	Insert Date ▼	Checked Out? ▼
Select Aliquot	AQID	Cage	Cane	Stack	Freezer	Original ID	Grid	Label Information	Volume (mL)	Insert Date	Checked Out?
<input type="checkbox"/>	001A-A01- <a href="#">Update</a>	1	01	A01	240B	00047303	A01	00047303	1.000	April 23, 2015, 5 p.m.	No
<input type="checkbox"/>	001A-A02- <a href="#">Update</a>	1	01	A01	240B	00047303	A02	00047303	1.000	April 23, 2015, 5 p.m.	No

- All fields of an Aliquot can be updated on this page.
- Required Fields have an \* after the label and if left blank a widget will notify the user and the submission will not go through.

Update an Aliquot

The current aliquot selected is 001A-A01

Lotid	Trialkey	Studyid	Studydesc	AliqId*	Originalid
20400047	Unknown	Unknown	ANA	001A-A01	0000047303
Solid*	Grid	Insertdate	Invnum		
001A	A01	2020-04-09 00:00:00	59		
Checkout*	Volume	Age	Gender		
1	,000	Unknown	Unknown		
Demogroup	Riskcategory	Riskgroupdesc	Bloodtype		
Unknown	Disease	ANA	Unknown		
Hiv1	Hiv2	Hiv3	Hbs	Hbc	Hav
Unknown	Unknown	Unknown	Unknown	Unknown	Unknown
Hcv	Mmr	Influenza	Sars cov2	Other	
Unknown	Unknown	Unknown	Unknown	ANA	
Manufacturer	Panelname	Paneldesc	Membernum		
Unknown	NA	NA	NA		
Samptypedesc	Samptypedesc	Dup	Isconfirmed		
Serum	SST	0	1		
Invcomment	Labelinfo	Notes	Reason*		
	00047303				

Submit

- Note that the Aliquot ID(Aqid) cannot be updated.
- If the user entered an incorrect Aqid then they must add a new aliquot with the correct Aqid and dispose of the incorrect Aqid.